President Grab Bars & Rails

Fixing and Maintenance Instructions

INTRODUCTION

Thank you for deciding to purchase from the Aidapt President Grab Rail/Grab Bar range. This range is manufactured from the finest materials and components available. When installed and used correctly they are designed to give many years of reliable trouble-free service. Please examine your product for any visible damage before using. If you see any damage or suspect a fault, please contact your supplier.

NB. This equipment must be installed by a competent person.

VY445/VY446/VY447	President White Grab Bar
VY448/VY449/VY450	President Blue Grab Bar
VY451/VY452/VY453	President Red Grab Bar
VY432/VY433	President White/Blue Grab Bar
VY425/VY426/VY427	President Ribbed Grab Bar (32 mm)
VY425B/VY426B/VY427B	President Ribbed Grab Bar (35 mm)

N.B. Any weight limits quoted are based on the product performance with secure fixings.

BEFORE USE

- Carefully remove all packaging. Avoid using any knives or other sharp instruments as this may damage the surface of the product.
- Check the product for any visible damage. If you see any damage or suspect a fault, do not use your product as it may compromise your safety, but contact your supplier for support.

FIXING GRAB BARS/GRAB RAILS

Firstly check the substrate where you intend to fix your Grab Rails/Grab Bars, ensuring that it is sound. Also check the type of substrate you intend to fix to, e.g. brick/ block timber or concrete.

Never drill through pre-stressed concrete to install a President product.

Always check the type of fixing you are going to use on this substrate for their compatibility.

We suggest that for sound brick and block that a minimum screw length of 50 mm x size 10 is used with a suitable rawl plug and screws, in accordance with the guidelines issued by the fixings manufacturer.

It is recommended that stainless steel screws are used so as to avoid the chance of rusting.

PLASTER BOARD (STUD TYPE) LATHE AND PLASTER

When attempting to fix to plasterboard stud or lathe and plaster, we suggest that you check the position of the studs and span them with a suitable timber plate of a minimum thickness of 40 mm x 150 mm high, securing the plate to the studs. After ensuring the soundness of the studs (using a suitable counter sunk screw), fix your grab bar to the timber plate.

Never install any grab bar or grab rail directly onto plaster board, MDF (medium density fibreboard), chipboard or other composite material. This will seriously compromise the safety of the user.

Always use all screw holes provided in the grab bar.

ELECTRICITY

Always check for concealed electrical cables when installing hinged arm supports. Also check your legal requirements under the NICE IC 16th Edition relating to Earth bonding of extraneous metal parts.

CLEANING

President Grab Rails/Grab Bars should be cleaned using a non-abrasive cleaner or mild detergent with a soft cloth. Only use a mild detergent or disinfectant; never use an abrasive cleaner. Do not exceed 50°C, and never autoclave Grab Rails/Grab Bars. Always ensure that you wipe dry all equipment and treat joints with a water dispersing agent, e.g. WD40.

REISSUE

If you reissue or are about to reissue this product, please thoroughly check all components for their safety.

If in any doubt, please do not issue or use, but immediately contact your supplier for service support.

CARE & MAINTENANCE

Please make a safety check of the product at regular intervals or if you have any concerns.

IMPORTANT INFORMATION

The information given in this instruction booklet must not be taken as forming part of or establishing any contractual or other commitment by Aidapt Bathrooms Limited, Aidapt (Wales) Ltd or its agents or its subsidiaries and no warranty or representation concerning the information is given.

Please exercise common sense and do not take any unnecessary risks when using this product; as the user you must accept liability for safety when using the product.

Please do not hesitate to contact the person who issued this product to you or the manufacturer (detailed below) if you have any queries concerning the assembly/use of your product.

Aidapt Bathrooms Ltd, Lancots Lane, Sutton Oak, St Helens, WA9 3EX

Telephone: +44 (0) 1744 745 020 • Fax: +44 (0) 1744 745 001 • Web: www.aidapt.co.uk

Email: accounts@aidapt.co.uk • adaptations@aidapt.co.uk • sales@aidapt.co.uk

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**Altai Germany (Altai Europe Ltd.)