This user manual contains essential information in relation to safe use and should be retained with the bed rail for reference

## Solo Bed Rail


Figure 1

**VY435**

## Fixing and Maintenance Instructions

EN

## INTRODUCTION

Thank you for deciding to purchase the Aidapt Solo Bed Rail. The Bed Rail is manufactured from the finest material and components available, when installed and used correctly it is designed to give many years of reliable trouble-free service. The Solo Bed Rail can be fitted to single, double, queen and king size beds. It is available in a tough epoxy coated finish. The solo bed rail is not designed to prevent users falling from their bed.

**⚠️ This equipment must be installed by a competent person and a risk assessment undertaken as to the suitability of the Bed Lever for the particular user and user's bed. Assessors should read and understand the MHRA publication; Bedrails: management and safe use, which can be downloaded from [www.gov.uk](http://www.gov.uk).**

**⚠️ 127 kg (20 st.) weight limit.**


**Do not exceed the weight limit stated – doing so could put the user at risk.**

## ⚠️ BEFORE USE

- Carefully remove all packaging. Avoid using any knives or other sharp instruments as this may damage the surface of the product.
- Check the product for any visible damage. If you see any damage or suspect a fault, do not use your product as it may compromise your safety, but contact your supplier for support.
- The Bed Transfer Aid hand grips have been designed to meet the latest MHRA guidelines on entrapment, however you must ensure that there is no risk to the user or carer after the installation and service of the Bed Transfer Aid.
- Do not allow children or unauthorised persons to play with the equipment or use the Bed Rail without proper supervision.

## ASSEMBLY AND INSTALLATION

1. Your Bed Rail is shipped ready assembled.
2. To install, remove the mattress from the bed where the Bed Rail is to be installed and check for any obstructions.
3. Place the Bed Rail onto the base of the bed with the large flat base on the bed frame and the looped handle facing upwards level with the side of the bed (see figures 1 and 2).
4. Move the Bed Rail up and down the bed to find the best position for the user to safely lift themselves.
5. Ensure that the frame of the Bed Rail is flush and square to the side of the bed.
6. Replace the mattress ensuring that the Bed Rail handle is pointing upwards and the unit is in the optimum position to suit the client's needs.

**Figure 2**

## ADDITIONAL SAFETY


For extra safety when in use, the Bed Rail can be fitted in conjunction with the strap and buckle set (available separately see figure 3).

To fit the strap and buckle set, please follow the instructions below:

1. Remove the mattress from the bed.
2. Turn the Bed Rail upside down so that the top of the handle and the end of the base are on the bed.
3. Pass the short strap (with the female buckle connector) under to the centre of the bar on the base of the Bed Rail closest to the handle and pass the buckle through the 'D' ring (see figure 4). Pull the buckle tight so that the strap is held securely in place against the frame.

**Figure 3**

**Figure 4**

4. Take the long strap and remove the mail buckle connector from the strap (see figure 5 and 6).

**Figure 5**

**Figure 6**

6. Feed the long strap under the centre of the bar on the base of the Bed Rail furthest from the handle and pass the plain end through the 'D' ring at the other end (see figure 7).
7. Feed the long strap under the centre of the bar on the base of the Bed Rail furthest from the handle and pass the plain end through the 'D' ring at the other end (see figure 7).
8. Re-attach the mail buckle connector to the long strap, turn the Bed Rail back up the right way and relocate it in its optimum position as set previously.


**Figure 7**

**Figure 8**

9. Pass the strap across the width of the bed frame, then underneath the bed and then up the other side.
10. Clip the two parts of the buckle together and tighten the long strap to remove any slack, ensuring that the strap is not twisted at any point (see figure 9).
11. Replace the mattress, check that the Bed Rail is secure and does not move by applying lateral pressure to the handle.

**Figure 9**

## CLEANING

The Bed Rail should be cleaned using a non-abrasive cleaner or mild detergent with a soft cloth only. Abrasive cleaners e.g. AJAX and/ or abrasive cleaning pads will seriously damage the product beyond repair and should not be used. Always ensure that you wipe dry the equipment after cleaning.

## REISSUE

If you reissue or are about to reissue this product, please thoroughly check all components for their safety.

This includes:

- Tightness of all nuts/bolts/castors and other screw in/bolt in/push in components.
- Upholstery (if applicable) for security, splits, exposed staples, etc.
- Straightness and integrity of all tubing and welds.

**If in any doubt, please do not issue or use, but immediately contact your supplier for service support.**

When reissuing, a risk assessment must be undertaken as to the suitability of the bed rail for the particular user and user's bed.

## CARE & MAINTENANCE

Please make a safety check of the product at regular intervals or if you have any concerns.

## IMPORTANT INFORMATION

The information given in this instruction booklet must not be taken as forming part of or establishing any contractual or other commitment by Aidapt Bathrooms Limited, Aidapt (Wales) Ltd or its agents or its subsidiaries and no warranty or representation concerning the information is given.

Please exercise common sense and do not take any unnecessary risks when using this product; as the user you must accept liability for safety when using the product.

Please do not hesitate to contact the person who issued this product to you or the manufacturer (detailed below) if you have any queries concerning the assembly/use of your product.

Manufactured by: Aidapt  
Bathrooms Limited Lancots Lane  
St Helens Merseyside, WA9 3EX  
UK

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**Altai Germany (Altai Europe Ltd.)**