



Wheelchair and Scooter Bags



**VA136ST/VA134ST/
VA136SS/VA132SS/VA134SS**

Fitting and Maintenance Instructions

INTRODUCTION

Thank you for deciding to purchase the Wheelchair or Scooter Bag from Aidapt. When used correctly they are designed to give many years of reliable, trouble-free service.

Please examine your product for any visible damage before using/installing. If you see any damage or suspect a fault, please do not use, but contact your supplier.

INTENDED USE

Wheelchair/ Scooter Bags are suitable for use with most standard wheelchairs/scooters. They provide the user and/or carer with extra storage.

FITTING INSTRUCTIONS

Wheelchair Bags

Simply hook a looped strap onto each of the handles to suspend the bag from the back of the wheelchair.

Scooter Bags

Simply hook the attached straps on the back of the Scooter Bag over the top of the seat with the bag on the back of the seat and adjust the horizontal strap as necessary for a snug and secure fit.

CLEANING

Clean your Wheelchair or Scooter Bag using a mild detergent with a slightly damp, soft cloth only. Only use a mild detergent or disinfectant; never use an abrasive cleaner.

After cleaning, leave the bag somewhere warm and dry to air. Do not leave in direct contact with a heat source.

CARE & MAINTENANCE

Please make a safety check of the product at regular intervals or if you have any concerns.

IMPORTANT INFORMATION

Please exercise common sense and do not take any unnecessary risks when using this product; as the user you must accept liability for safety when using the product.

Please do not hesitate to contact the person who issued this product to you or the manufacturer (detailed below) if you have any queries concerning the assembly/use of your product.

**Aidapt Bathrooms Ltd, Lancots Lane, Sutton Oak, St Helens, WA9
3EX**

**Telephone: +44 (0) 1744 745 020 • Fax: +44 (0) 1744 745 001 • Web:
www.aidapt.com**

Email: sales@aidapt.co.uk